



PL 94-104 Łódź, ul. Obywatelska 187/189  
tel. +48 42 640 48 70, fax +48 42 688 33 84



## Instructions for use of modeling wax

### Application:

- to use in making dentures.

### Properties:

- very good plasticity achieved by using beeswax and natural resins in its production.

### Direction for use:

Before use plastify the plate of wax by dipping it in warm water or holding over the flame of a burner. The plate of wax, prepared in above mentioned method can use in prosthodontic.

### Storage:

Keep packaged, in a dry place, protected from the sun and far from any source of heat.

**Best before:** 3 years from the date of production.

### Different types of modelling wax:

- HARD
- SOFT